
11/96-Pr-Sf-

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The Compton scattering equipment set makes possible qualitative and quantitative observation of the Compton effect. Compton scattering illustrates the particle character of γ radiation particularly well, and thus plays an important role as an introductory experiment to quantum mechanics and the dualism of wave and particle. The experiment setup allows you to measure the changes in the wavelength of the scattered γ radiation as a function of the scattering angle and intensity. This in turn permits determination of the differential effective cross-section.

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1 Equipment description, scope of delivery, technical data; (see Fig. 1)

- ① Sample holder, i.d. 12 mm, with collimator
- ② Detector holder with lead shielding
- ③ Moveable additional shielding
- ④ Aluminum scatterer
- ⑤ Experiment panel with angular scale and polar coordinates

- ① Sample holder as lead block with centered hole, 12 mm dia., for accommodating radioactive preparations. Additional blind hole for inserting one of the steel pins (included in scope of delivery) as directional indicator for collimator hole.
- ② Detector holder with lead shielding for defined direction of incoming γ radiation, designed for the scintillation counter (559 90) in conjunction with the detector output stage (559 91). Additional blind hole which for inserting one of the steel pins (included in scope of delivery) as directional indicator for detector window.
- ③ Moveable additional shielding: lead block serves to reduce the intensity of unscattered γ radiation, particularly for small scattering angles and short distances between source, scatterer and detector.
- ④ Aluminum scatterer, cylindrical pure aluminum rod as center of scattering.

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2 Operation

- First connect MCA CASSY (529 780) with the high-voltage power supply, 1.5 kV (521 68), scintillation counter (559 90) and the detector output stage (559 91) in accordance with the respective instruction sheets.
- Insert the scintillation counter into the detector holder ② from the rear.
- Set up the detector holder ② and the sample holder ① on a sheet of polar-coordinate paper (for easy angle determination) as shown in Fig. 2.
- Insert the ^{137}Cs preparation (559 809) (in the sample holder ①) from the rear.
- By varying the high voltage level on the high-voltage power supply (521 68), adjust the gain so that the photo peak of the ^{137}Cs can be recorded with no problem.
- Using the ^{137}Cs (559 809) and ^{241}Am preparations, record the calibration spectra at $\vartheta = 0$. Calibrate the energy with the program "MCA".
- For every desired scattering angle ϑ , record a spectrum for ^{137}Cs . Use additional shielding as necessary.
- For a reliable determination of the energy of the scattering peak E_ϑ , a measurement with and without scatterer should be made for every angle ϑ to permit evaluation of the difference of the spectra. The determination of this differential is essential for determining the peak intensity I_ϑ .
- Determine the energy of the scatter radiation E_ϑ , and thus the change in the wavelength, from the position of the peak.
- Comparison with theory, in accordance with Fig. 2:

$$E_\vartheta = E_0 \cdot \frac{1}{(1 + \alpha \cdot (1 - \cos \vartheta))}$$

where

$$\alpha = \frac{E_0}{m_e \cdot c^2}$$

E_0 : Energy of unscattered γ radiation
 $m_e c^2$: Rest mass of the electron

- The intensity I_ϑ of the scattered radiation is obtained through integration of the photo peak.
- The relative intensities I_ϑ of the scatter radiation peaks can be compared with the predictions of the Klein-Nishina formula for the differential effective cross section $d\sigma/d\Omega$. The calibration factor C can be calculated using the formula:

$$C = \frac{1}{n} \cdot \sum_{\vartheta=0}^n \frac{I_\vartheta}{\frac{d\sigma}{d\Omega}}$$

σ : cross section
 Ω : solid angle of leakage radiation
 n : number of different scattering angles.

Notes

- 1) The measuring time per scattering angle must be selected according to the desired accuracy and the distances between source, scatterer and detector.
- 2) Greater distances between source, scatterer and detector increase the angular resolution, but require longer measuring times due to the lower counting rates.
- 3) At small scattering angles, particularly for short distances, additional lead shielding is required to suppress the background caused by unscattered γ radiation.


Fig. 2


Fig. 3


Fig. 4.1


Fig. 4.2


Fig. 4.3


Fig. 4.4

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Differential spectra for Compton effect with ¹³⁷Cs (559 809) for various angles ϑ , determined from the spectra recorded with and without aluminium