## **FÍSICA EXPERIMENTAL IV**

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### Balmer Lamp Power Supply Unit for Balmer Lamp Balmer Lamp, Deuterated

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The Balmer lamp (451 13) provides - in conjunction with a spectrometer array of the type used in schools - the four visible lines  $H_\alpha$ ,  $H_\beta$ ,  $H_\gamma$  and  $H_\delta$  of the hydrogen spectrum (Balmer series). The quantitative evaluation of the spectrum permits determination of the wavelengths, so that the Balmer series formula

$$\frac{1}{\lambda} = R \left( \frac{1}{2^2} - \frac{1}{n^2} \right) \quad n = 3, 4, 5, 6$$

$$R = 3,29 \cdot 10^{15} \text{ s}^{-1}$$

(Rydberg frequency)

can be confirmed by experiment.

If a high-resolution spectrometer (resolution approx. 0.1 nm) is available, the line doublets of a hydrogen/deuterium mixture can be demonstrated using the Balmer lamp, deuterated (451 41).

Literature:

Physics Experiments (Cat. No. 599 861, in German): Atomic and Nuclear Physics, experiment 6.2.2.1

## 1 Safety Notes

- Always disconnect the mains plug from the wall outlet before making any changes to the Balmer experiment setup - particularly when changing the Balmer lamp!
- Always disconnect the mains plug from the wall outlet before changing the fuse as described in section 4.
- Caution: do not touch the Balmer lamp - very hot!

## 2 Description, Technical Data

### 2.1 Balmer lamp (451 13)

The Balmer lamp is an AC operated gas discharge tube filled with water vapor. The fused tube is supplied with water vapor by means of a water supply bound in hygroscopic material. The water molecules are broken down by the electric discharge into atomic hydrogen and a hydroxyl group. A high-temperature resistant capillary inside the lamp confines the discharge within a narrow space so that a high concentration of atomic hydrogen is realized. This atomic hydrogen is responsible for the intense Balmer spectrum; interference due to molecular hydrogen bands cannot occur.

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An oxidant, assisted by suitable catalysts, ensures that hydrogen formed during operation is oxidized to water so that the water inside the tube is recycled. This recycling of water involves the deposition of reddish-brown metal oxides which, however, will not have a negative effect on the capillary part of the tube.

n	Linie	Wellenlänge $\lambda = \frac{c}{f}$
3	H $\alpha$	656,28 nm (rot)
4	H $\beta$	486,13 nm (türkis)
5	H $\gamma$	434,05 nm (blau)
6	H $\delta$	410,17 nm (violett)

n	Line	Wavelength $\lambda = \frac{c}{f}$
3	H $\alpha$	656.28 nm (red)
4	H $\beta$	486.13 nm (turquoise)
5	H $\gamma$	434.05 nm (blue)
6	H $\delta$	410.17 nm (violet)

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Table 1  
Wavelengths of the lines observable with the Balmer lamp

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## 2.2 Balmer lamp, deuterated (451 41)

In the Balmer lamp, deuterated, the water filling consists of about 10 % deuterated water. The individual lines of the Balmer spectra of hydrogen and deuterium differ slightly in position.

Linien	Linienabstand $\Delta\lambda$ (nm)
H $\alpha$ - D $\alpha$	0,177
H $\beta$ - D $\beta$	0,131
H $\gamma$ - D $\gamma$	0,117
H $\delta$ - D $\delta$	0,111

Line	Line spacing $\Delta\lambda$ (nm)
H $\alpha$ - D $\alpha$	0.177
H $\beta$ - D $\beta$	0.131
H $\gamma$ - D $\gamma$	0.117
H $\delta$ - D $\delta$	0.111

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Table 2  
Line spacing of doublets in the spectrum of the Balmer lamp, deuterated

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\*) We would like to thank Dr. Klein of the Physical Institute of the University of Cologne for his kind permission to reproduce the spectra shown here. These were recorded using a high-resolution Heath scanning monochromator, model EU 700, with RCA photomultiplier 931 A.

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### 2.3 Power supply unit for the Balmer lamp (451 14)

The power supply (a) F2 consists of the power supply unit operated from the AC mains outlet and the Balmer lamp socket, which is permanently connected with this unit, contact-protected and high voltage-proof. The high AC voltage required to operate the Balmer lamp is produced in the high-voltage power supply by means of a high-reactance transformer.

The high-voltage cable from the rear of the power supply unit ends at a distributor box (h) (see fig. 3). The tubes (g) fitted to the distributor box are used for two purposes:

- 1) to support the Balmer lamp, and
- 2) to enclose the high-voltage connections.

At the ends of the tubes (g) are the plastic lamp sockets. The lower socket (d) is permanently connected with the tube (g) while the upper socket (e) can be demounted by detaching the clip (f). The lamp holder can be attached via the stand rod (c) either on the portable high-voltage power supply (clamp with knurled screw (b)) or to the rider of an optical bench or to stand material.

The on-off switch with operating indicator lamp and the mains-voltage selector switch (k) with integrated fuse holder (see Fig. 5) are located on the rear of the device.

#### Technical data

Mains connection voltage: 230 V AC; 50/60 Hz, set at works;  
convertible to 115 V AC  
No-load voltage: approx. 3500 V<sub>rms</sub>  
Operating data with lamp connected: arc voltage approx. 1500 V  
arc current approx. 50 mA  
Fuses: for 230 V AC: T 1.25 B (fitted at works)  
for 115 V AC: T 2.5 D (included in scope of supply)


Fig. 3

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#### 3.1 Preparing the Balmer lamp (see fig. 3)

Important!

When setting up the Balmer lamp, make sure the power supply is not connected to the mains.

1. Fit the lamp holder with stand rod ③ to clamp ② of the high-voltage power supply.
2. Detach the clip ① at the upper lamp socket ⑤, pull it off downward and place it on the distributor box ④.
3. Remove the top lamp socket ⑤ by pulling it upward. The high-voltage cable with contact spring bushing is exposed.
4. Insert the aluminium socket of the Balmer lamp into the center of the lower, permanently mounted lamp socket ④ and press it down. Make sure that the metal pin of the lamp is seated in the contact spring bushing of the socket.
5. Connect the cable, led from the tube ⑥ via the contact spring bushing.
6. Press on the upper lamp socket ⑤ until part of the ribbing in tube ⑥ is covered.
7. Fix the upper lamp socket ⑤ to tube ⑥ by means of clip ①.

#### 3.2 Putting the Balmer lamp into operation

1. Set the high-voltage power supply to the correct mains voltage by means of the selector switch.
2. Connect to the mains.
3. Set the switch on the high-voltage power supply to "ON". The ON state is indicated by a red lamp above the switch. The Balmer lamp is ignited at once.

#### 3.3 Notes on operation

1. To ensure steady burning of the Balmer lamp, a certain operating temperature is required (approx. 55 °C, measured at the aluminium socket). If the operating temperature is too low (e.g. immediately upon switching on), discharges outside the capillary may occur, causing the lamp to flicker. Steady discharge conditions will be reached after an operating time of approx. 10 to 15 minutes.

If at high ambient temperatures the operating temperature becomes too high (above approx. 70 °C, measured at the aluminium socket), the lamp may go out after a longer period of operation. In this case the lamp can be normally operated again after it has cooled down.

2. If after a long period of operation the hot Balmer lamp does not ignite after being switched off and immediately switched on again, wait until it has cooled down before switching it on once more.
3. The lamp should only be operated in the vertical position, with the aluminium socket facing downward. In other operating positions the light will start to flicker after some minutes or will even go out.
4. At initial operation, a band-like background will appear in addition to the strong Balmer spectrum. This background will disappear after approx. one hour's operation. The oxygen lines in the red and green spectral region, which are very weak compared to the Balmer lines, will remain (cf. Fig. 1).
5. With hot Balmer lamps the yellow sodium line will occasionally be visible in the middle section of the capillary. In this case it is advisable to observe the upper and lower ends of the capillary where the disturbing yellow line will not appear.

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#### 3.4.1 Spectrum of the Balmer lamp (451 13)

Subjective observation:

Pocket spectroscope	467 02
or	
School spectroscope	467 112
or	
Spectroscope and goniometer	467 23
or	
Copy of a Rowland grating, approx. 6000 lines/cm	471 23

Objective observation (Projection arrangement see Fig. 4):

1 Copy of a Rowland grating	471 23
1 Optical bench e.g	460 43
1 Bench clamp	301 06
1 Adjustable slit	460 14
1 Holder with spring clips	460 22
1 Translucent screen	441 53
1 Lens in frame, $f = 50$ mm	460 02
1 Lens in frame, $f = 100$ mm	460 03
6 Leybold multiclamps	301 01

Set up the experiment as shown in Fig. 4, but initially without the grating.  
In a darkened room, focus the Balmer lamp sharply on the slit using the lens  $f = 50$  mm. Focus the light from the slit sharply on the translucent screen by moving the lens  $f = 100$  mm back and forth; then attach the grating to its holder using the spring clips.  
Open the slit until a reasonable compromise is reached between brightness and separation (resolution) of the lines.


Fig. 4

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### 3.5 Spectrum of the Balmer lamp, deuterated (451 41)

Minimum resolution of spectral apparatus for separation of the line doublets: 0.1 nm (e.g. Heath scanning monochromator, model EU 700 and photomultiplier RCA 931 A).

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
#### 4 Replacing the fuse (Fig. 5)

Important:

- Disconnect the mains plug.
- Check the fuse for the correct value: T 1.25 B for 230 V AC  
T 2.5 D for 115 V AC

To exchange the primary fuse, (required when a fuse blows or when converting to 115 V AC) put a coin into slot ① of the voltage selector and fuse holder and turn it until the white mark ② is next to position "0" as shown in Fig. 5. Catch the fuse in your hand as it is pushed out of opening ③ by the spring. Insert the new fuse into opening ③ and press down by means of a pointed object (e.g. ball pen, screw-driver), at the same time turning and pressing down the coin in the slot ①.

Set the voltage so that the respective value of 220 V (for AC mains voltage of 230 V) or 110 V (for AC mains voltage of 115 V) is opposite the white mark ②.

Fig. 5

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### 1 Description

The copy of a Rowland grating is copied from an original Rowland grating on a transparent collodium foil and is about 23 mm x 35 mm in size. The grating has about 600 lines per mm; the grating constant, the distance between two adjacent lines, is therefore 1/600 mm. The foil is inserted in 50 mm x 50 mm metal frame between two glass plates.

The collodium foil is subject to minor changes of shape due to meteorological influences. Therefore, the spectroscopical resolution is not quite the same as in case of the original grating.

### 2 Use

The frame is clamped into the holder with spring clips (460 22) and either mounted on a saddle base (300 11) or on an optical bench (e.g. 460 43).


### 3 Experiments

a) Projection of diffraction spectra  
In order to produce brilliant diffraction spectra for group viewing we recommend using the mercury vapour lamp (451 15) in an experimental arrangement of the type shown in Fig. 2. The lens,  $f = 200$  mm, forms an image of the slit on a wall 2 m to 3 m away. After mounting the grating directly behind the projecting lens the spectral lines of mercury can be seen on the wall in a bright display.

Fig. 2

- ① 451 15/19
- ⑤ 460 22
- ⑩ 460 04
- ⑪ 460 43
- ⑬ 460 14
- ⑭ 471 23
- ⑰ 460 02
- ⑲ 451 30

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b) Individual viewing of spectra using spectrum lamps

The spectra of spectrum lamps (451 011 - 451 111 with 451 16 and universal choke 451 30, see relevant Instruction Sheet) are best viewed in an arrangement of the type shown in Fig. 3. If the eye is brought up close to the grating, spectra of several orders can be viewed towards both sides.


Fig. 3

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c) Measurements of wavelengths or of grating constants

These measurements are e.g. done in an experimental arrangement as shown in Fig. 4, with a spectral tube which supplies a line spectrum with defined wave lengths. Please see to it that the grating is aligned in parallel to the rule at a distance of about 70 cm to 120 cm from it. The spectral tube must be directly in front of the rule.

The angle of diffraction follows from the distance between the rule and the grating and from the deviation of the spectral lines that can be read off on the rule. If the wavelengths are known, the grating constant must be computed and if the grating constant is known, the wavelengths can be determined.



Fig. 4

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## Spectrum Lamps Ne, He, Cd, Hg, Hg-Cd, Tl, Na Housing for Spectrum Lamps

Fig. 1

The spectrum lamps are AC-operated gas-discharge lamps emitting the line spectrum of metal vapors or noble gases respectively with high spectral purity and high luminous intensity.

The spectrum lamps are set up in the housing for spectrum lamps (451 16); the power lead of the housing is connected to the universal choke (451 30), which supplies precisely the current needed to operate the lamps.

#### Examples of Experiments:

- Determination of wavelengths
- Inversion of the Na-D line (using 451 111)
- Investigation of UV radiation (using 451 062)
- Photoelectric effect (using 451 062)
- Determination of  $h$  (using 451 062)

## 1 Safety notes

Spectral lamps (451 011 and following numbers) in housing (451 16) should be connected to the mains only via the universal choke (451 30)

Keep a minimum distance of 3 cm between light exit opening and optical element (e.g. diaphragm, lens) to avoid overheating due to accumulation of heat.

## 2 Technical Data

### 2.1 Spectrum lamps (451 011 to 111)

- ① Glass bulb (10 cm x 2 cm dia.) with burner
- ② 9-pin pico socket matched to housing (451 16)

Gas filling*)	Ne	He	Cd	Hg	Hg-Cd	Tl	Na
Cat. No.	451 011	451 031	451 041	451 062	451 071	451 081	451 111
Operating voltage	30 V	60 V	15 V	45 V	30 V	15 V	15 V

Operating current:  $1 \text{ A} \pm 10\%$

(ensured with supply via universal choke, 451 30, at mains voltages between 200 V AC and 230 V AC)

\*) Spectrum lamps filled with caesium, potassium, rubidium and zinc are also available to special order.

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## 2.2 Housing for spectrum lamps (451 16)

- ③ Lower part of housing with 9-pin pico socket, built-in starter
- ④ Upper part of housing, movable on ③, with light emission aperture (2.5 cm dia.)  
Dimensions of housing: 16 cm x 7 cm dia.
- ⑤ Stand rod, 5.5 cm long, 10 cm dia.
- ⑥ Connecting lead, 80 cm long, with multipin plug, to be connected, safe to touch, to the universal choke (451 30) by means of union nut.

## 3 Operation

*Additionally required:*

Universal choke 451 30  
Optical elements and  
Stand material, e.g. according to Fig. 2

Plug in spectrum lamp, connect housing with universal choke, switch on.

Allow to warm up, the full radiation output is achieved only after some minutes, except spectrum lamp He (451 031) and spectrum lamp Ne (451 011).

In case of power failure, allow the warm spectrum lamp Hg (451 062) to cool down some minutes before switching it on again (otherwise it will not ignite).


Fig. 2